



DETERMINATION OF BUPRENORPHINE IN HUMAN PLASMA BY XLC-MS/MS USING SYMBIOSIS™ PHARMA

December 2005
0053.063-01

APPLICATION INFO

Introduction

Symbiosis™ Pharma is Spark Holland's unique solution for integrated online SPE-LC-MS automation (XLC-MS). The system offers large flexibility in processing different types of samples selecting one of the three fully automated operational modes LC-MS; XLC-MS; AMD (Advanced Method Development).

This application will demonstrate the compatibility of the Symbiosis™ Pharma to add internal standard during a XLC-MS routine with excellent accuracy, precision and linearity over the whole calibration range.

Buprenorphine (BUP) is a powerful semi-synthetic analgesic. BUP's molecular structure is related to Morphine, but it is estimated to be at least 25 times more potent. Initially, it was available for the treatment of acute and chronic pain. (Brand names are: Temgesic®, Buprenex® and Subutex®).




Figure 1: Buprenorphine, CAS# 52485-79-7, C₂₉H₄₁NO₄, Mw 467.6472. Physical properties: Water solubility 0.63 mg/L, LogP 4.98, pKa dissociation constant 8.31.

Method Development



Figure 2: Symbiosis™ Pharma System

The AMD mode of Symbiosis Pharma in conjunction with the HySphere method development cartridge tray (Spark p.n. 0722.650) enables "quick sorbent

screening" for most suitable SPE cartridge and optimal wash conditions for clean-up. The following data was obtained in less than 1 hour using generic pre-defined SPE conditions of the Symbiosis™ Pharma.




Figure 3: Chromatogram of Buprenorphine in human plasma during sorbent screening.

After the sorbent screening the C18HD cartridge is selected and the wash procedure is optimized.




Figure 4: Method Development Cartridge Tray

After optimizing the wash protocol an overall recovery of >95% was calculated from the peak area of a neat solution and serum sample (figure 4)




Figure 5: Overlay of 10 ng/mL plasma XLC sample and 10 ng/mL aqueous LC sample.

For more information visit our website:
[HTTP://WWW.SPARKHOLLAND.COM](http://www.sparkholland.com)

Spark
HOLLAND

XLC-MS Protocol

Autosampler Conditions

30 µL sample is injected using µL pick up mode injection routine. In this injection mode the autosampler adds the internal standard (Buprenorphine D4) to sample.




Figure 6: Flow path of Symbiosis Pharma with I.S. addition

Washing is performed with two wash solvents;
 Wash solvent 1: 50% ACN with 0.1% Formic Acid.
 Wash solvent 2: 90% ACN.

Wash solvent	Wash volume	Valve wash
1	700 µL	no
2	700 µL	no
1	700 µL	yes
2	700 µL	yes
1	1400 µL	yes

Table 1: Autosampler wash routine.

SPE conditions

Cartridge:	10 x 2 mm HySphere C18HD (Spark PN:0722.609)	
Solvation:	1 mL ACN	5 mL/min
Equilibration:	1 mL 5% ACN in 1% NH ₄ OH pH 9.5	5 mL/min
Sample Loading:	0.5 mL 40% ACN in 1% NH ₄ OH pH 9.5	1 mL/min
Washing:	1 mL 40% ACN in 1% NH ₄ OH pH 9.5	5 mL/min
Elution	2 min. with LC Gradient	
Total SPE time is 2 min 30 sec.		

Table 2: SPE conditions.

LC conditions

Column:	Waters Xterra® MS C18 2.1 X 50 mm 3.5µ
Mobile phase A:	0.1% Formic Acid in Water
Mobile phase B:	0.1% Formic acid in ACN

Time (mm:ss)	Flow (mL/min.)	A (%)	B (%)
00:01	0.30	90	10
00:05	0.30	90	10
02:35	0.30	10	90
02:50	0.30	10	90
03:00	0.30	90	10
05:00	0.30	90	10




Figure 7: LC gradient

MS Conditions

A Sciex API 3000 with a Turbo IonSpray in positive mode is used.

	Buprenorphine	Buprenorphine D4 (IS)
Q1 mass	468.45	475.45
Q3 mass	55.10	59.10
Dwell time	150	150
DP	86	86
FP	270	270
CE	103	103
CXP	4	4

Table 3: Compound dependable MS settings. (MS parameters: NEB = 15; CUR = 13; IS = 5500; TEM = 400; CAD = 6; EP = 10)

Result

The following samples are prepared in Human Plasma. The autosampler adds the internal standard during the injection routine (Vial-to-File™ concept).

- Calibration standards: 0.1; 0.5; 1.0; 5.0; 10; 50; 100 ng/mL
- QC samples: 0.1; 1; 10; 100 ng/mL

Chromatograms

Figure 8 and 10 are representing chromatograms of the upper and lower limits of the calibration curve. The carry-over represented in the chromatogram of the blank (figure 9) is calculated as <0.02% (< 20% of LOQ).




Figure 8: Chromatogram representing 100 ng/mL Buprenorphine.




Figure 9: Chromatogram representing a blank after 100 ng/mL.




Figure 10: Chromatogram representing 0.1 ng/mL Buprenorphine.

Linearity, Accuracy and Precision

A calibration curve was determined by combining the results of 6 repeated injections of a full set of calibration standards. This resulted in a R^2 of 0.997 with 1/X Weighting.




Figure 11: Peak area vs. concentration from Buprenorphine (corrected for Internal Standard) from the 6 combined set of calibration standards.

Sample (ng/mL)	CV (%)	Accuracy (%)
0.1	4.78	105
1.0	1.92	93.4
10	1.13	89.1
100	2.84	100

Table 4: Accuracy and Precision calculated from 4 QC series

Reproducibility

In this application, the internal standard is added by the autosampler.

To determine the reproducibility, the area of the internal standard of 105 samples are displayed in figure 12. The calculated %RSD of the internal standard is 3.9%




Figure 12: Peak Area of internal standard in plasma RSD = 3.9%

Conclusions

From this study it is concluded that within a time frame of two days it is possible to develop a XLC-MS method with an absolute recovery of >95%, and run 6 sets of calibration standards with a linear range from 0.1-100 ng/mL ($R^2 = 0.998$ with a 1/X weighting). The developed method has an accuracy between 85-105% and a precision of <5% CV.

Furthermore it shows that the Reliance™ autosampler in Symbiosis™ Pharma is able to add internal standards with an excellent reproducibility of 3.9% RSD.

About Spark

Since 1982 Spark has provided the HPLC and LC/MS markets with state-of-the-art autosamplers, column ovens and sample preparation solutions. Solid Phase Extraction with on-line elution into HPLC and LC/MS systems was pioneered by Spark and introduced in the early 90's. Spark, ISO 9001 certified, does basic research, product development, production, sales and marketing in-house, guaranteeing quality from start to finish. With 25% of the employees working in research and development Spark continues to invest in the future, making sure we can deliver the solutions you need to improve your business results. Innovation and quality are keywords when talking about our development efforts.

Spark System Solutions BV
Bendienplein 5
7815 SM Emmen, the Netherlands

P +31 591631700
F +31 491645900
E Solutions@Sparkholland.com
W www.sparkholland.com