


## On-line versus off-line solid phase extraction for the determination of Microcystins in Water by LC-MS/MS

1. Spark Holland B.V. Bendienplein 5, Emmen, the Netherlands 2. AGBAR Aigües de Barcelona. Av. Diagonal 211 Barcelona, Spain

### OVERVIEW

- Microcystins (MC) are natural toxins produced by Cyanobacteria in natural waters.
- Its toxicity has been demonstrated in both animals and human beings. The can produce extensive liver damage and promote liver tumor formation in subacute exposures.
- The World Health Organization (WHO) suggested a maximum daily intake concentration of MC of 1µg/L for adults. The European Water Framework Directive (2000/60/EC) has established a maximum permitted level of MC in drinking water at the same level suggested by the WHO, 1µg/L

| Compound | X          | Y       | CAS Nr.      | Nominal Mass |
|----------|------------|---------|--------------|--------------|
| MC-LR    | Leucin     | Arginin | 101043-37-2  | 994          |
| MC-YR    | Tryptophan | Arginin | 1010-64-48-6 | 1044         |
| MC-RR    | Arginin    | Arginin | 111755-37-4  | 1037         |



### INTRODUCTION

The present work describes the first fully automated method, based on on-line Solid Phase Extraction (SPE) – ultra fast liquid chromatography (UHPLC)-tandem mass spectrometry (MS/MS) for the determination of microcystins (MC-LR, -YR, -RR) in surface waters. Combination of both on-line SPE and UHPLC has been achieved by using special fused core™ chromatographic columns which allowed maintaining the high efficiencies obtained with UHPLC reducing the system backpressure. The analysis of MC was then performed in a total analysis time of 15 minutes. High recoveries, ranging from 99% to 115% were obtained and the focusing mode was used during SPE elution in order to focus analytes in the column head and to avoid peak broadening and high efficiencies were then achieved (N=34000-80000). The quality parameters of the method were established; limits of quantification (LOQ) were similar than those obtained with the off-line SPE ranging from 6 to 10 ng/L while lower run-to-run and day-to-day precisions were obtained in the range of 1.8-2.9% and 2.3-3.6% respectively.

### Objectives

- Develop an on-line SPE-LC method coupled to tandem mass spectrometry for the determination of MC in water resources.
- Maintain or improve the benefits achieved with the existing off-line SPE UPLC-MS/MS in terms of resolution and sensitivity.
- Improve the method in terms of cost- and time efficiency.
- Perform the XUHPLC-MS/MS coupling and validate the method.

### EXPERIMENTAL CONDITIONS

#### ONLINE SPE-LC (SYMBIOSIS™ PICO, SPARK HOLLAND)

##### AUTOSAMPLER

Injection volume: 10 mL  
 Sample: Water (+5% Methanol pH=7)  
 Needle wash: 700 µL, 40% ACN, 0.2% FA

##### LC

Column: Ascentis Express C18 (Supelco, Sigma)  
 Column oven: 50 °C  
 Mobile phase: A) ACN, 0.1% FA B) Water, 0.1% FA

##### Gradient 1:

| time (m:s) | flow (µL/min) | % A | % B |
|------------|---------------|-----|-----|
| 0.01       | 400           | 5   | 95  |
| 2.00       | 400           | 5   | 95  |
| 2.02       | 500           | 25  | 75  |
| 7.45       | 500           | 95  | 5   |
| 8.15       | 500           | 95  | 5   |
| 8.30       | 500           | 5   | 95  |
| 10.00      | 500           | 5   | 95  |




### SPE extraction

Cartridge: HySphere C18-HD cartridge (10 x 2mm ID-7µm, Spark)

Conditioning solvents:  
 1 mL MeOH (5mL/min)  
 2 mL Water (5mL/min)


Sample loading: with Water  
 Washing:  
 1 mL MeOH:Water (20:80)

Focusing (Elution):  
 HPD syringe: MeOH 200 µL 100 µL/min  
 Gradient pump: Mobile phase (A:B; 5:95) 400µL/min



### MS : AB/Sciex API 3200 Qtrap MS/MS

|                  | RT (min) | Q1   | DP(V) | Q3  | CE(V) | Q3  | CE(V) |
|------------------|----------|------|-------|-----|-------|-----|-------|
| MC-RR            | 3.46     | 995  | 116   | 135 | 87    | 105 | 127   |
| MC-YR            | 3.85     | 1045 | 146   | 135 | 99    | 105 | 129   |
| MC-LR            | 3.97     | 1038 | 131   | 135 | 103   | 127 | 117   |
| Nodularin (I.S.) | 3.68     | 825  | 96    | 103 | 127   |     |       |




### Method Performance

|       | LOD (ng/L) | LOQ (ng/L) | Run-to-run (%RSD)* | Day-to-day (%RSD)* | Linearity |
|-------|------------|------------|--------------------|--------------------|-----------|
| MC-RR | 1.2        | 10         | 2.9                | 3.6                | 0.995     |
| MC-YR | 0.4        | 6          | 1.8                | 2.3                | 0.997     |
| MC-LR | 0.8        | 8          | 2.2                | 2.3                | 0.998     |

- Satisfactory LODs were obtained in the low ng/L level.
- Excellent linearity's were achieved for each compound.

### Off-line SPE VS On-line SPE

|                       | Off-line       | On-line    |
|-----------------------|----------------|------------|
| Sample volume         | 100 mL         | 10 mL      |
| SPE Solvent volume    | 12 mL          | 1.2 mL     |
| Evaporation           | N <sub>2</sub> | -          |
| Reconstitution volume | 6 mL           | -          |
| Total Analysis time   | 2 hours        | 15 minutes |



### Off-line SPE VS On-line SPE

|       | SPE-(Off-line) UPLC- MS/MS |     |       |       | SPE-(On-line) LC-MS/MS |     |       |       |
|-------|----------------------------|-----|-------|-------|------------------------|-----|-------|-------|
|       | Rt (min)                   | Rs  | N     | Asym. | Rt (min)               | Rs  | N     | Asym. |
| MC-RR | 1.83                       | -   | 35130 | 0.8   | 3.46                   | -   | 33965 | 1.1   |
| MC-YR | 2.35                       | 4.0 | 21114 | 1.0   | 3.85                   | 1.8 | 79985 | 1.3   |
| MC-LR | 2.46                       | 0.6 | 23260 | 1.1   | 3.97                   | 0.8 | 52755 | 1.1   |

### DISCUSSION and CONCLUSION

- The benefits of both on-line SPE system and UHPLC technologies (coupled to tandem MS) have been successfully combined for the analysis of microcystins in waters.
- The combination of both techniques was performed by using a fused core particle column which allowed to obtain high efficiencies and reduced backpressures.
- The three MC were separated in less than 4 minutes without losing resolution.
- Focusing mode during the SPE elution step avoided peak broadening and the high efficiencies achieved with the UHPLC were maintained
- Under these conditions the total analysis time was reduced to 15 minutes yielding to a reduction in both time and solvent consumption without losing efficiency.
- High recoveries (99-115%) and low LOQs (6-10 ng/L) were obtained for all the compounds.

### Off-line SPE routine:

For the off-line SPE method, a Waters Acquity ultra-performance liquid chromatograph (UPLC®) system, equipped with a quaternary pump system (Milford, MA, USA) using an Acquity BEH C18 column (5.0 cm x 2.1 mm i.d., 1.7 µm particle size) [Waters Corp., Milford, MA, USA] was used.

Waters Oasis® HLB Cartridges were rinsed with 10 mL of methanol and 10 mL of Milli-Q water, samples were percolated through the cartridge at a flow rate of 10 mL/min. The SPE cartridge was washed with 8 mL of a 20% methanol aqueous solution and dried with nitrogen gas for 10 minutes. Analytes were eluted from the cartridge using 6 mL of methanol and the extracts were stored at -20 °C. Prior to injection, the extract was evaporated to dryness at 40 °C under a stream of nitrogen in a TurboVap LV evaporator from Zymark (Zymark, Hopkinton, MA). Samples were reconstituted in 200 µL of an aqueous solution with a 5% of methanol. The final extracts were filtered through 0.2 µm before injection

### AgBar, Barcelona, Spain

The Agbar Group is at the head of a big holding made up of more than 150 companies and 140 years of history which operates in all the fields related to community services: complete water cycle and health, among others. The Agbar Group is a world leader both in the fields which constitute the core of its activities related to the water cycle and in the business lines demanded by the new society, particularly in the services of health. Its international expansion, with presence in United Kingdom, Chile, China, Algeria, Colombia, Cuba and Mexico consolidates its leadership and its potential of growth as a group rendering service to more than 37 million inhabitants all over the world. With a workforce over 13,000 employees, operating revenues of €2,861 million and consolidated income attributed to the parent company of €352.5 million in 2007, the Agbar Group has established continuous rates of growth and of value creation which allow it to advance day by day, adapting to the needs of a constantly evolving world, and also to occupy a leading position in most of the areas in which it develops its activities.

