

A new automated approach for the determination of peptides in human plasma using on-line SPE – LC – MS/MS

Spark Holland, Pieter de Keyserstraat 8, 7800AJ, Emmen, The Netherlands, Phone (+31) 591 631 700, E-mail: emile.koster@sparkholland.com, Website: www.sparkholland.com

Overview

- A new automated sample extraction technique for the determination of peptides through Mixed Mode Strong Anion exchange is developed.
- The new method is tested with 3 peptides with different molecular masses: 2100, 5900 and 6200 amu.

Introduction

Sample preparation for analytes in biological matrices is a critical part of successful LC MS/MS analyses. Analysis of peptides is particularly challenging due to some unique properties such as size, carry-over, polarity differences and poor solubility, which can affect extraction performance. All these properties have a severe impact on the recovery of large peptides from biological matrices like plasma and serum. A new extraction approach for peptides is evaluated: integrated clean-up and analysis of peptides using an on-line SPE anion exchange procedure with a mixed mode strong anion exchange cartridge. Several peptides drugs from different suppliers were investigated with this procedure.

Experimental conditions

Samples were prepared by spiking pooled human plasma with a peptide. Prior to extraction, all plasma samples were diluted 1:4 with water. 20 µL Methanol sample (100 µL) was directly placed into the autosampler and analyzed fully automated. Online SPE is performed on a disposable cartridges; each sample uses a fresh SPE cartridge. After extraction, the peptides are eluted from the cartridge onto a reverse phase analytical column utilizing the HPD focus mode.

Online –SPE-LC (Symbiosis™ Pico, Spark Holland)

LC

Column: Phenomenex Jupiter 4µ Proteo 50x2 mm [p/n.00B-4396-00]
 Mobile phase 1: 0.1% Formic Acid in water
 Mobile phase 2: 95% Methanol and 5% 0.1% Formic Acid in water


SPE

Cartridge: Oasis® MAX 30 µm 10x1 mm (Waters corporation)
 Conditioning: 1 mL Methanol at 5 mL/min
 Equilibration: 1 mL 20% Methanol in 1% NH4OH water at 5 mL/min

Sample injection: 100 µL spiked plasma
 Sample extraction: 200 µL Methanol in 1% NH4OH water at 0.5 mL/min

Cartridge wash 1: 1 mL 20% Methanol in 1% NH4OH water at 5 mL/min
 Cartridge wash 2: 1 mL 80% Methanol in 1% NH4OH water at 5 mL/min

HPD Focus elution: 300 µL Methanol at 100µL/min
 Clamp flush: 1 mL 80/20 acetonitrile/water 0.2% FA at 5 mL/min



ESI-MS (API4000, Applied Biosystems)

ESI-MS/MS conditions (positive mode)

MS-settings					
Compound Mass	Q1 mass	Q3 mass	DP	CE	CXP
2100	498.5	110.2	76	101	20
5900	1184.2	454.4	51	57	14
6200	1239.7	731.6	86	51	22

General settings: IS 500 TEM 400 CAD7 CUR 15 GS1 75 OS2 40 EP 10 Dwell 100

Results and Discussion

The sample clean-up is optimized by increasing the percentage organic modifier in the wash solvent. Figure 1 shows the wash optimization for the 6200 AMU peptide. The other two used peptides showed similar results.




Figure 1: wash optimization for peptide 6200. peptide 5900 and 2100 showed same results

As shown in figure 1 the mixed mode anion exchange cartridge allows for a high organic wash without reducing the analyte recovery. For optimal clean-up a first wash step with 20% ACN was programmed to remove interfering proteins (no denaturing on SPE cartridge), followed by a second high organic wash at 80% ACN to remove non-polar matrix compounds (see chromatograms below).




Figure 2: 6200 AMU peptide. Spiked plasma sample extracted on Oasis MAX (10X1) mm




Figure 3: 5900 AMU peptide. Spiked plasma sample extracted on Oasis MAX (10X1) mm




Figure 4: 2100 AMU peptide. Spiked plasma sample extracted on Oasis MAX (10x1) mm

As shown in figure 2-4 clean and reproducible chromatograms were obtained. Subsequently, calibration curves were determined by injecting a full set of calibration standards ranging from 0.5 to 500 ng/mL spiked human plasma standards. The calibration curves are calculated with a 1/X weighting. The curves for the peptides 6200, 5900 and 2100 are displayed below (r > 0.996 for all three peptides).




Figure 5: 6200 AMU peptide. Spiked plasma calibration curve from 0.5 to 500 ng/mL




Figure 6: 5900 AMU peptide. Spiked plasma calibration curve from 0.5 to 500 ng/mL




Figure 7: 2100 AMU peptide. Spiked plasma calibration curve from 0.5 to 500 ng/mL

The method was evaluated using spiked human plasma samples. The developed XLC-MS/MS method has an absolute recovery of more than 90%. The accuracy is ranging from 89-117% and the precision 2.1-7.9 %CV. The estimated quantitation limits for the peptides are more than sufficient to allow the analytical method to be used for therapeutic drug monitoring.

Conclusion

- A simple sample preparation method was developed for the determination of the peptides in plasma.
- In this concept 100 µL of a spiked plasma sample is loaded on a mixed mode anion exchange cartridge.
- The developed concept of automated on-line SPE – LC – MS/MS delivers clean and reproducible results.
- Good linearity was obtained for all 3 peptides in human plasma over the 0.5-500 ng/mL range (r > 0.994).
- More than 100 µL sample could be loaded onto the cartridge to achieve even lower detection limits.
- The Symbiosis™ Pico method is currently used routinely in several pharmaceutical labs (various companies) for analyzing different size of peptides (< 10K amu) with similar results.

